



# Laboratory Report

UL  
85 John Road  
Canton, MA 02021 USA  
Phone: (781)821-2200  
Fax:(781)8219266  
Report#: 1002845163  
PAGE: 1 of 2

GREG BARROW  
LiveDo USA, Inc.  
4925 LIVEDO DR WILSON NC 27893 USA

Received Date: Sep 12, 2025  
Test Date: Sep 12, 2025 - Sep 19, 2025  
Report Date: Sep 22, 2025

Item Description	LivDry Overnight Protective Underwear	Supplier	LiveDo USA Inc.
Client Reference	0032320	Portal Id	GCTN013415
Lot #	U2 250709 2373 EXP 07/28	Inventory Id	8912407
UPC	892684001226	Order No	16012492

**Approved By**  
DANIEL PROIA  
Client Services Manager II

LETTERS & REPORTS: UL Verification Services, Inc. (UL) letters and reports are issued for the exclusive use of the clients to whom they are addressed. No quotations from reports or use of the UL name is permitted except as expressly authorized in writing. Letters and reports apply only to the specific materials, products or processes tested, examined or surveyed as received and are not necessarily indicative of the qualities of apparently identical or similar materials, products or processes. The liability of UL with respect to services rendered shall be limited to the amount of consideration paid for such service and not include any consequential damages. This report or certificate does not relieve sellers/suppliers from their contractual responsibility with regard to the quality/quantity of this delivery, nor does it prejudice clients' right to claim towards sellers/suppliers for compensation for any apparent and/or hidden defects not detected during our random inspection or testing. UL has not performed a complete analysis of the product. The results contained in this report indicate that the product has passed or failed the specific tests only, as determined by client requirements. These test results, even if rated as "Passed," do not indicate or certify that the product is safe for commercial or consumer use.



# Laboratory Report

UL  
 85 John Road  
 Canton, MA 02021 USA  
 Phone: (781)821-2200  
 Fax:(781)8219266  
 Report#: 1002845163  
 PAGE: 2 of 2

Criteria	LivDry Overnight Protective Underwear
<b>Method: Total Absorbency and Retention Capacity of Absorbent Products [ISO 11948-1 and NWSP 070.10.R1]</b>	
Dry Weight, Replicate 1	93.93 g
Dry Weight, Replicate 2	93.67 g
Dry Weight, Replicate 3	93.05 g
Dry Weight, Replicate 4	92.43 g
Dry Weight, Replicate 5	93.14 g
Dry Weight, Average	93.24 g
Dry Weight, Minimum	92.43 g
Dry Weight, Maximum	93.93 g
Total Absorbency, Replicate 1	2063.36 g
Total Absorbency, Replicate 2	2074.25 g
Total Absorbency, Replicate 3	2103.88 g
Total Absorbency, Replicate 4	2048.69 g
Total Absorbency, Replicate 5	2100.44 g
Total Absorbency, Average	2078.12 g
Total Absorbency, Minimum	2048.69 g
Total Absorbency, Maximum	2103.88 g

LETTERS & REPORTS: UL Verification Services, Inc. (UL) letters and reports are issued for the exclusive use of the clients to whom they are addressed. No quotations from reports or use of the UL name is permitted except as expressly authorized in writing. Letters and reports apply only to the specific materials, products or processes tested, examined or surveyed as received and are not necessarily indicative of the qualities of apparently identical or similar materials, products or processes. The liability of UL with respect to services rendered shall be limited to the amount of consideration paid for such service and not include any consequential damages. This report or certificate does not relieve sellers/suppliers from their contractual responsibility with regard to the quality/quantity of this delivery, nor does it prejudice clients' right to claim towards sellers/suppliers for compensation for any apparent and/or hidden defects not detected during our random inspection or testing. UL has not performed a complete analysis of the product. The results contained in this report indicate that the product has passed or failed the specific tests only, as determined by client requirements. These test results, even if rated as "Passed," do not indicate or certify that the product is safe for commercial or consumer use.